

| | | |
|---|---------------------|----------------------------|
| 1 | Status | Letter Ballot |
| 2 | Date of Last Update | 2024/11/11 |
| 3 | Person Assigned | David Clunie |
| 4 | | mailto:dclunie@dclunie.com |
| 5 | Submitter Name | David Clunie |
| 6 | | mailto:dclunie@dclunie.com |
| 7 | Submission Date | 2024/05/18 |

| | | |
|----|---|--|
| 8 | Correction Number CP-2415 | |
| 9 | Log Summary: Add PAXgene fixation system | |
| 10 | Name of Standard | |
| 11 | PS3.16 | |
| 12 | Rationale for Correction: | |
| 13 | PAXGene is a proprietary fixation reagent used in preparation of pathology specimens. | |
| 14 | Correction Wording: | |

Amend DICOM PS3.16 as follows (changes to existing text are bold and underlined for additions and ~~struckthrough~~ for removals):

CID 8114 Specimen Fixative

Version: ~~20190124~~yyyymmdd

Table CID 8114. Specimen Fixative

| Coding Scheme Designator | Code Value | Code Meaning | SNOMED-RT ID | UMLS Concept Unique ID |
|--------------------------|-----------------------|-------------------------------------|--------------|------------------------|
| SCT | 434162003 | Neutral Buffered Formalin | C-2141C | C0492002 |
| SCT | 433474002 | Bouin's fluid | F-62235 | C0053963 |
| SCT | 431510009 | Formalin | C-2141B | C0949307 |
| SCT | 433338005 | Carnoy's fluid | F-62231 | C2317379 |
| SCT | 434295000 | Formol sublimate | F-62238 | C0621539 |
| SCT | 433471005 | Helly's fluid | F-62233 | C2317380 |
| SCT | 430028007 | Michel's medium | F-6220F | C1550080 |
| SCT | 433473008 | Zenker's fluid | F-62234 | C2317478 |
| SCT | 52836003 | Paraformaldehyde | C-21403 | C0070066 |
| SCT | 2869004 | Acetic acid | C-21624 | C0000983 |
| SCT | 259153006 | Chloroform | C-20830 | C0008238 |
| SCT | 430821002 | Chromium trioxide | C-12916 | C0055630 |
| SCT | 419442005 | Ethanol | C-21047 | C0001962 |
| SCT | 111095003 | Formaldehyde | C-21402 | C0016564 |
| SCT | 11496005 | Mercuric chloride | C-13321 | C0025417 |
| SCT | 259221006 | Methanol | C-2102B | C0001963 |
| SCT | 13931001 | Osmium tetroxide | C-15211 | C0029385 |
| SCT | 24215009 | Picric acid | C-21919 | C0071044 |
| SCT | 19893005 | Potassium dichromate | C-13518 | C0032829 |
| <u>NCIt</u> | <u>C185113</u> | <u>PAXgene Tissue System</u> | | <u>C5669715</u> |